

PRODUCT BROCHURE

M.E.D.
SURGICAL





M.E.D. Surgical, CELEBRATING 25 YEARS OF HEALTHCARE EXCELLENCE.

Established in 1991, M.E.D. Surgical specialised in the distribution of OR products. The company has grown significantly since then and our mission is to provide best in-class service across all divisions of expertise. Our focus is to ensure all our customers receive a level of service that not only meets but exceeds their expectations.

An extensive portfolio covering a range of specialties including Ophthalmology, Infection Prevention, Anaesthetics, Cardiology/Vascular, Urology/Gynaecology and Rehabilitation, has made us the leading provider of specialist medical devices in Ireland.

Our commitment to excellent customer service has nurtured and developed multiple long-standing partnerships across public and private hospitals and is exemplified by the employment of 33 people nationwide, including 21 dedicated sales specialists.

M.E.D. Surgical provide the highest level of efficiency, turnaround times and close working partnerships with their customers. A policy of continuous self-appraisal and attention to detail ensure the ongoing expansion of our customer base.



“M.E.D. Surgical provide the highest level of efficiency, turnaround times & close working partnerships with their customers.”

ORTHOPAEDIC SOFT GOODS AND BRACING



M.E.D. Surgical are delighted to have announced their latest addition to their portfolio. The Orthopaedic Soft Goods and Bracing division was launched in May 2017 and primarily focuses on A&E, Plaster Room, Physiotherapy, Rehab Units, Orthopaedic Out patients, Minor Injuries Unit and Tissue Viability. The Orthopaedic Soft Goods and Bracing Division is set to provide for all Hospitals, The Private Sector, G.P's and Patients Direct in the Republic of Ireland.

Products available are all upper and lower limb Orthoses. Covering all Trauma Departments, Rehabilitation Units, Orthotics, and Post-Surgical Units. This division represents a number of leading manufacturers of Orthopaedic Soft Goods and Bracing such as:

- **Talarmade**
Specialists in Orthopaedic Soft Goods
Based in the UK
- **Danmic**
Specialists in Hand Therapy and
Rehabilitation Products
Based in USA

M.E.D. Surgical also provide a dedicated sales team to drive this division. With over 30 years' experience in Soft Goods and Post Op Bracing between UK and Ireland, we offer our customers on-going support and assistance on a daily basis. In addition to this we also provide in house Training, Seminars and Educational Days working alongside Clinicians and Company Specialists.

M.E.D. Surgical offer a full service to all hospitals enabling Purchasing Departments to rationalise their suppliers using M.E.D. Surgical as a one stop shop, therefore saving resource on procurement.

TALARMADE

DELUX ARM SLING

FEATURES & BENEFITS

- Breathable lightweight polyester cotton minimises skin irritation
- Soft padded plush foam straps for comfort with easy touch-and-close fastenings
- Thumb loop and full circumferential waist strap to provide stabilisation for the arm

RECOMMENDED FOR

- Immobilisation of forearm, elbow and arm
- Post upper limb trauma
- Post upper limb surgery



SIZE	INCHES	CMS	PRODUCT CODE
SMALL	14	36	DAS/S
MEDIUM	16	41	DAS/M
LARGE	19	48	DAS/L
X-LARGE	21	52	DAS/XL

MULTI ARM SLING

FEATURES & BENEFITS

- Brushed velour with fully bound edges
- Fully adjustable shoulder strap and full circumference torso strap
- Universal sizing

RECOMMENDED FOR

- Immobilisation of forearm, elbow and arm



SIZE	INCHES	CMS	PRODUCT CODE
UNIVERSAL	-	-	PAS/U

VERSATILE ARM SLING

FEATURES & BENEFITS

- High quality foam sling
- This product is designed for single patient use only
- Washing or laundering is not recommended
- Only for use in contact with intact skin
- Available in 2 sizes

RECOMMENDED FOR

- Pre or post surgery where the forearm needs to be held in an elevated position



SIZE	UNIT OF SALE	PRODUCT CODE
SMALL / CHILD	BOX OF 10	VAS/C
ADULT	BOX OF 10	VAS/A

HAND POSITIONING BRACE

FEATURES & BENEFITS

- Malleable frame allows static progressive splinting
- No tools or heat required to adjust
- Padded frame has removable washable cover
- Cover lined with synthetic sheepskin and velour

RECOMMENDED FOR

- Post surgery / Injury positioning
- Wrist drop and finger flexor tightness
- Contractures following stroke, arthritis & trauma
- Not recommended when moderate to severe spasticity is present



WRIST CIRCUMFERENCE			LEFT	RIGHT
SIZE	INCHES	CMS	PRODUCT CODE	PRODUCT CODE
SMALL	5.5 - 8	14 - 20.5	HPB/SL	HPB/SR
MEDIUM	6.75 - 9.25	17 - 23.5	HPB/ML	HPB/MR
LARGE	8 - 10.5	20.5 - 26.5	HPB/LL	HPB/LR



PALM PROTECTOR

FEATURES & BENEFITS

- Synthetic sheepskin to improve comfort
- Closed cell foam separator
- Touch and close fastening
- Available in 1 size (left or right)

RECOMMENDED FOR

- Severe MCP and finger flexion contractures
- Improving hand hygiene



		LEFT	RIGHT
SIZE	UNIT OF SALE	PRODUCT CODE	PRODUCT CODE
ONE SIZE	SINGLE	BMPP/L	BMPP/R

ELASTIC WRIST BRACE

FEATURES & BENEFITS

- Woven elastic material for comfort and durability with bound edges providing compressive support
- Removable anatomically formed aluminium palmar bar individually mouldable for custom fit
- Simple touch-and-close fastening with circumferenti support strap for wrist
- Available in single units and packs of 6

RECOMMENDED FOR

- Mild wrist sprain and strains
- Mild symptoms of Carpel Tunnel syndrome
- Rheumatoid / Osteoarthritis
- Healing Fracture



WRIST CIRCUMFERENCE			LEFT	RIGHT
SIZE	INCHES	CMS	PRODUCT CODE	PRODUCT CODE
SMALL	5.5 - 6.75	14 - 17	EWB/SL	EWB/SR
MEDIUM	6.25 - 7	16 - 18	EWB/ML	EWB/MR
LARGE	7 - 7.75	18 - 20	EWB/LL	EWB/LR
X-LARGE	7.5 - 8.25	19 - 21	EWB/XLL	EWB/XLR

For a pack of 6, add '6' at the end of the product code

AIRPRENE WRIST BRACE

FEATURES & BENEFITS

- High grade breathable Airprene material with stretch fabric on both sides with bound edges providing compressive support to wrist
- Removable anatomically formed aluminium palmar bar individually mouldable for custom fit
- Simple touch-and-close fastenings with circumferential support strap for wrist
- Available in single units and packs of 6

RECOMMENDED FOR

- Mild wrist sprain and strains
- Mild symptoms of Carpel Tunnel syndrome
- Rheumatoid / Osteoarthritis
- Healing fracture



WRIST CIRCUMFERENCE			LEFT	RIGHT
SIZE	INCHES	CMS	PRODUCT CODE	PRODUCT CODE
X-SMALL	4.25 - 5.5	12 - 14	AWB/XSL	AWB/XSR
SMALL	5.5 - 6.75	14 - 17	AWB/SL	AWB/SR
MEDIUM	6.25 - 7	16 - 18	AWB/ML	AWB/MR
LARGE	7 - 7.75	18 - 20	AWB/LL	AWB/LR
X-LARGE	7.5 - 8.25	19 - 21	AWB/XLL	AWB/XLR

For a pack of 6, add '6' at the end of the product code

BASIC WRIST BRACE

FEATURES & BENEFITS

- Removable aluminium palmar bar easily shaped to wrist
- Simple touch-and-close fastening
- Available in single units and packs of 6

RECOMMENDED FOR

- Mild wrist sprains
- Mild symptoms of carpal tunnel syndrome
- Rheumatoid / Osteoarthritis



WRIST CIRCUMFERENCE			LEFT	RIGHT
SIZE	INCHES	CMS	PRODUCT CODE	PRODUCT CODE
X-SMALL	4.25 - 5.5	12 - 14	AWB/XSL	AWB/XSR
SMALL	5.5 - 6.75	14 - 17	AWB/SL	AWB/SR
MEDIUM	6.25 - 7	16 - 18	AWB/ML	AWB/MR
LARGE	7 - 7.75	18 - 20	AWB/LL	AWB/LR
X-LARGE	7.5 - 8.25	19 - 21	AWB/XLL	AWB/XLR

For a pack of 6, add '6' at the end of the product code

COOL COMFORT WRIST BRACE

FEATURES & BENEFITS

- Pre-bent palmar bar allows 15° wrist extension
- Easy to fit and remove

RECOMMENDED FOR

- Osteoarthritis
- Rheumatoid Arthritis
- Sprains
- Strains
- Carpal Tunnel Syndrome



WRIST CIRCUMFERENCE			LEFT	RIGHT
SIZE	INCHES	CMS	PRODUCT CODE	PRODUCT CODE
X-SMALL	UP TO 6	UP TO 15	CCWB/XSL	CCWB/XSR
SMALL	6 - 6.5	15 - 17	CCWB/SL	CCWB/SR
MEDIUM	6.5 - 7.5	17 - 19	CCWB/ML	CCWB/MR
LARGE	7.5 - 8.5	19 - 22	CCWB/LL	CCWB/LR
X-LARGE	8.5 +	22+	CCWB/XLL	CCWB/XLR

WRIST ULNAR DEVIATION BRACE

FEATURES & BENEFITS

- A uniquely functioning wrist brace for the treatment of UP deformity
- The breathable neoprene materials used along with the towelling lining, allows for reduction in heat build and patient perspiration
- Contains a malleable aluminium palmar stay, individually moulded for custom fit

RECOMMENDED FOR

- Early MCP ulnar deviation deformity
- Rheumatoid Arthritis



WRIST CIRCUMFERENCE			LEFT	RIGHT
SIZE	INCHES	CMS	PRODUCT CODE	PRODUCT CODE
X-SMALL	4.5 - 5	11 - 13	WUDB/XSL	WUDB/XSR
SMALL	5.5 - 6	14 - 15	WUDB/SL	WUDB/SR
MEDIUM	6.5 - 7	16 - 18	WUDB/ML	WUDB/MR
LARGE	8 - 8.5	19 - 22	WUDB/LL	WUDB/LR
X-LARGE	9 - 9.5	23 - 24	WUDB/XLL	WUDB/XLR

DELTA ELASTIC WRIST BRACE

FEATURES & BENEFITS

- Seam direction improves patient comfort and fit
- High quality aluminium stay
- Available in 2 lengths, 17cms and 21cms
- Available in single units and packs of 6

RECOMMENDED FOR

- Mild wrist sprain and strains
- Mild symptoms of Carpal Tunnel Syndrome
- Rheumatoid Osteoarthritis
- Healing fracture



DELTA 17CMS STANDARD LENGTH

WRIST CIRCUMFERENCE			LEFT	RIGHT
SIZE	INCHES	CMS	PRODUCT CODE	PRODUCT CODE
X-SMALL	4.5 - 5	11 - 13	DWB/XSL	DWB/XSR
SMALL	5.5 - 6	14 - 15	DWB/SL	DWB/SR
MEDIUM	6.5 - 7	16 - 18	DWB/ML	DWB/MR
LARGE	8 - 8.5	19 - 22	DWB/LL	DWB/LR
X-LARGE	9 - 9.5	23 - 24	DWB/XLL	DWB/XLR

For a pack of 6, add '6' at the end of the product code

DELTA PLUS 21CMS LONG LENGTH

WRIST CIRCUMFERENCE			LEFT	RIGHT
SIZE	INCHES	CMS	PRODUCT CODE	PRODUCT CODE
X-SMALL	4.5 - 5	11 - 13	DPWB/XSL	DPWB/XSR
SMALL	5.5 - 6	14 - 15	DPWB/SL	DPWB/SR
MEDIUM	6.5 - 7	16 - 18	DPWB/ML	DPWB/MR
LARGE	8 - 8.5	19 - 22	DPWB/LL	DPWB/LR
X-LARGE	9 - 9.5	23 - 24	DPWB/XLL	DPWB/XLR

For a pack of 6, add '6' at the end of the product code

AIRPRENE WRIST THUMB BRACE

FEATURES & BENEFITS

- High grade breathable Airprene material with stretch fabric on both sides with bound edges providing compressive support to wrist & thumb
- Removable anatomically formed aluminium palmar bar and thumb stay individually mouldable for custom fit
- Simple touch-and-close fastenings
- Available in single units and packs of 6

RECOMMENDED FOR

- Mild wrist sprain and strains
- Mild symptoms of Carpal Tunnel Syndrome
- Rheumatoid / Osteoarthritis
- Healing Fracture



WRIST CIRCUMFERENCE			LEFT	RIGHT
SIZE	INCHES	CMS	PRODUCT CODE	PRODUCT CODE
SMALL	5.5 - 6.75	14 - 17	AWTB/SL	AWTB/SR
MEDIUM	6.25 - 7	16 - 18	AWTB/ML	AWTB/MR
LARGE	7 - 7.75	18 - 20	AWTB/LL	AWTB/LR
X-LARGE	7.5 - 8.25	19 - 20	AWTB/XLL	AWTB/XLR

For a pack of 6, add '6' at the end of the product code

VARIABLE COMPRESSION WRIST/THUMB SPICA

FEATURES & BENEFITS

- Provides variable compression to the wrist
- Adjustable allowing the user to tighten or loosen the support
- Reinforced metal stay provides additional stabilisation
- Elastic strap provides additional compression

RECOMMENDED FOR

- Rheumatoid Arthritis
- Carpal Tunnel Syndrome
- Tenosynovitis



SIZE	INCHES	CMS	PRODUCT CODE
UNIVERSAL	-	-	BMVCWTS/U

ELASTIC WRIST THUMB BRACE

FEATURES & BENEFITS

- Woven elastic material for comfort and durability with bound edges providing compressive support
- Removable anatomically formed aluminium palmar bar and thumb stay individually mouldable for custom fit
- Simple touch-and-close fastenings
- Available in single units and packs of 6 (Add 6 at the end of the product code)

RECOMMENDED FOR

- Mild wrist / thumb sprain and strains
- Stabilising hyper or hypo tonic wrist / thumb flexors
- Rheumatoid / Osteoarthritis
- Healing Fracture



WRIST CIRCUMFERENCE			LEFT	RIGHT
SIZE	INCHES	CMS	PRODUCT CODE	PRODUCT CODE
SMALL	5.5 - 6.75	14 - 17	EWTB/SL	EWTB/SR
MEDIUM	6.25 - 7	16 - 18	EWTB/ML	EWTB/MR
LARGE	7 - 7.75	18 - 20	EWTB/LL	EWTB/LR
X-LARGE	7.5 - 8.25	19 - 21	EWTB/XLL	EWTB/XLR

For a pack of 6, add '6' at the end of the product code

AIRPRENE THUMB SPICA

FEATURES & BENEFITS

- High grade breathable Airprene material with stretch fabric on both sides and bound edges providing compressive support to wrist & thumb
- Metal stay extends along thumb to restrict movement at the CMC and MCP joints leaving the IP joint free
- Universal left or right fitting
- Available in single units and packs of 6 (Add 6 at the end of the product code)

RECOMMENDED FOR

- Thumb injury requiring immobilisation
- CMC ligament sprains & strains



SIZE	INCHES	CMS	PRODUCT CODE
SMALL	5 - 6	13 - 15	ATS/S
MEDIUM	6 - 6.75	15 - 17	ATS/M
LARGE	6.75 - 7.5	17 - 19	ATS/L
X-LARGE	7.5 - 8.25	19 - 21	ATS/XL

For a pack of 6, add '6' at the end of the product code

ELASTIC THUMB SPICA

FEATURES & BENEFITS

- Metal stay extends along thumb to restrict movement at the CMC and MCP joints leaving the IP joint free
- Available in single units and packs of 6

RECOMMENDED FOR

- Thumb injury requiring immobilisation
- CMC ligament sprains & strains



WRIST CIRCUMFERENCE			LEFT	RIGHT
SIZE	INCHES	CMS	PRODUCT CODE	PRODUCT CODE
SMALL	5 - 6	12.5 - 15	ETS/SL	ETS/SR
MEDIUM	6 - 6.75	15 - 17	ETS/ML	ETS/MR
LARGE	6.75 - 7.5	17 - 19	ETS/LL	ETS/LR
X-LARGE	7.5 - 8.75	19 - 22.5	ETS/XLL	ETS/XLR

For a pack of 6, add '6' at the end of the product code

COOL COMFORT THUMB ABDUCTION SPLINT

FEATURES & BENEFITS

- Wrap-around support made from breathable neoprene with terry cloth lining for patient comfort
- Includes firm mouldable insert inside thumb strap
- Thumb splint helps to abduct and stabilise the first metacarpal and open the thumb web space

RECOMMENDED FOR

- Thumb deformity and contracture due to neurological conditions
- Median nerve injuries
- Arthritis



WRIST CIRCUMFERENCE			LEFT	RIGHT
SIZE	INCHES	CMS	PRODUCT CODE	PRODUCT CODE
SMALL	6 - 7	15 - 18	CCTAS/SL	CCTAS/SR
MEDIUM	7 - 8	18 - 20	CCTAS/ML	CCTAS/MR
LARGE	8 - 9	20 - 23	CCTAS/LL	CCTAS/LR

COOL COMFORT THUMB RESTRICTION SPLINT

FEATURES & BENEFITS

- Wrap-around support made from breathable neoprene with terry cloth lining for patient comfort
- Provides gentle compression to relieve discomfort associated with CMC Arthritis or RSI
- Neoprene strap wraps around the thumb CMC joint and the first metacarpal to direct CMC joint support

RECOMMENDED FOR

- Osteoarthritis
- Medial subluxation of the thumb
- CMC joint



WRIST CIRCUMFERENCE			LEFT		RIGHT	
SIZE	INCHES	CMS	PRODUCT CODE	NHSSC CODE	PRODUCT CODE	NHSSC CODE
SMALL	6 - 7	15 - 18	CCTRS/SL	GQT1697	CCTRS/SR	GQT1698
MEDIUM	7 - 8	18 - 20	CCTRS/ML	GQT1695	CCTRS/MR	GQT1696
LARGE	8 - 9	20 - 23	CCTRS/LL	GQT1693	CCTRS/LR	GQT1694

POLYTHENE THUMB ORTHOSIS

FEATURES & BENEFITS

- Plastic shell holds thumb in exact degree of flexion
- Elastic strap provides comfortable fit during movement
- Easily trimmed or heated for exact fit

RECOMMENDED FOR

- DeQuervain's Tenosynovitis
- Gamekeeper's thumb / UCL tear / skiers thumb
- Stabilisation of the 1st MP joint



WRIST CIRCUMFERENCE			LEFT	RIGHT
SIZE	INCHES	CMS	PRODUCT CODE	PRODUCT CODE
SMALL	5 - 6	13 - 15	PTO/SL	PTO/SR
MEDIUM	6 - 7	15 - 18	PTO/ML	PTO/MR
LARGE	7 - 8	18 - 20	PTO/LL	PTO/LR

DIGITEC MALLET FINGER SPLINT

FEATURES & BENEFITS

- Made of high quality rigid plastic
- Supports DIP joint in a neutral position
- 5 piece set

RECOMMENDED FOR

- Immobilisation of the DIP joint



SIZE	CIRCUMFERENCE OF DIP (CMS)	PRODUCT CODE
1	5	DTMFS/1
2	5.5	DTMFS/2
3	5.8	DTMFS/3
4	6.1	DTMFS/4
5	6.3	DTMFS/5
5.5	6.8	DTMFS/5.5
6	7.3	DTMFS/6
7	8.4	DTMFS/7

DIGITEC FINGER SPLINT

FEATURES & BENEFITS

- Nylon shell
- EVA foam lined

RECOMMENDED FOR

- Protection of distal finger
- Immobilisation of acute injury



LENGTH OF SPLINT		
SIZE	CMS	PRODUCT CODE
SMALL	5.6	DTFS/S
MEDIUM	6.6	DTFS/M
LARGE	7.6	DTFS/L
X-LARGE	9.2	DTFS/XL

3 PANEL MATERNITY SUPPORT

FEATURES & BENEFITS

- Provides variable tension via the use of individual straps
- Firm circumferential support around the pelvis holding the two sacroiliac joints and symphysis pubis joint

RECOMMENDED FOR

- Relief from diagnosed ante natal or post natal back pain



WAIST CIRCUMFERENCE			
SIZE	INCHES	CMS	PRODUCT CODE
XX-SMALL	35 - 38	90 - 96	3PMS/1
X-SMALL	38 - 40	96 - 102	3PMS/2
SMALL	40 - 42	102 - 108	3PMS/3
MEDIUM	42 - 46	108 - 116	3PMS/4
LARGE	46 - 49	116 - 126	3PMS/5
X-LARGE	50 - 53	126 - 136	3PMS/6
XX-LARGE	53 - 57	136 - 146	3PMS/7
XXX-LARGE	57 - 61	146 - 156	3PMS/8

UNIVERSAL MATERNITY BELT

FEATURES & BENEFITS

- Low profile elastic perinatal sacroiliac support
- Easy touch and close fastening
- One size fits all
- Posterior stays for comfort preventing roll-over

RECOMMENDED FOR

- Low back pain associated with pregnancy



SIZE	INCHES	CMS	PRODUCT CODE
UNIVERSAL	-	-	MB-W/U

ECONOMY 3 PANEL KNEE IMMOBILISER

FEATURES & BENEFITS

- Versatile 3 panel design
- Made from high quality foam laminate with towelling lining
- Fully bound edges for durability
- Hook closure allows easy positioning of the detachable side panels allowing any size of leg to be accommodated
- Touch-and-close through loop fasteners for a secure fit
- 4 aluminium stays may be contoured for custom fit and all aluminium stays fully removable for easy customisation and to aid laundering
- Available in 9 lengths

RECOMMENDED FOR

- Patella dislocations
- Ligament injuries
- Pre operative pain conditions

LENGTH			
SIZE	INCHES	CMS	PRODUCT CODE
12	12	30	EC-3PKI/12
14	14	35	EC-3PKI/14
16	16	40	EC-3PKI/16
18	18	45	EC-3PKI/18
20	20	51	EC-3PKI/20
22	22	56	EC-3PKI/22
24	24	61	EC-3PKI/24
26	26	66	EC-3PKI/26
28	28	71	EC-3PKI/28



ANKLEGUARD STIRRUP ANKLE BRACE WITH AIR/GEL PADS

FEATURES & BENEFITS

- Removable double welded Air/Gel bladders
- May be heated or cooled to provide hot or cold therapy
- Streamlined design aids fit in shoes for early protected weight-bearing
- Universal stirrup shape fits left or right

RECOMMENDED FOR

- Stabilisation of chronic ankle instability
- Moderate to severe sprains
- Malleolar fractures
- Conservative and post-operative care following ligament injury / surgery

SIZE	INCHES	CMS	PRODUCT CODE
UNIVERSAL	-	-	AGAGSB/U



AIRPRENE KNEE BRACE

FEATURES & BENEFITS

- Front opening wrap-around design
- Made from high quality 4 way stretch breathable neoprene laminate with soft feel nylon fabric
- Fully bound edges
- Popliteal cut out to back of knee for added patient comfort
- Removable dual pivot polycentric hinge with hyperextension stop
- Posterior strap adjustment allows for one time adjustment and proper positioning of hinges
- Touch-and-close fastening with easy-grip tabs
- Universal fits left or right leg
- Available in beige or black

RECOMMENDED FOR

- Osteoarthritis
- Moderate ligament instabilities
- Traumatic or degenerative deformities of knee joint
- Unstable valgus or varus knee



AIRPRENE KNEE BRACE - BEIGE

KNEE CIRCUMFERENCE			SHORT 30CM	LONG 40CM
SIZE	INCHES	CMS	PRODUCT CODE	PRODUCT CODE
SMALL	11.75 - 13.5	30 - 34	ASKB/S	AVKB/S
MEDIUM	13.5 - 15	34 - 38	ASKB/M	AVKB/M
LARGE	15 - 16.5	38 - 42	ASKB/L	AVKB/L
X-LARGE	16.5 - 18	42 - 46	ASKB/XL	AVKB/XL

AIRPRENE KNEE BRACE - BLACK

KNEE CIRCUMFERENCE			SHORT 30CM	LONG 40CM
SIZE	INCHES	CMS	PRODUCT CODE	PRODUCT CODE
SMALL	11.75 - 13.5	30 - 34	ASKB-B/S	AVKB-B/S
MEDIUM	13.5 - 15	34 - 38	ASKB-B/M	AVKB-B/M
LARGE	15 - 16.5	38 - 42	ASKB-B/L	AVKB-B/L
X-LARGE	16.5 - 18	42 - 46	ASKB-B/XL	AVKB-B/XL

CANVAS POST OP SHOE

FEATURES & BENEFITS

- Durable soft padded upper and tongue
- Return through loop touch-and-close fasteners accommodates the oedematous or bandaged foot comfortably
- Lightweight flexible EVA sole with ribbed non-slip finish
- Universal fits right or left foot

RECOMMENDED FOR

- Post surgery requiring accommodation of bandages
- Oedematous foot



SIZE	SOLE LENGTH CMS	UK SHOE SIZE	PRODUCT CODE
X-SMALL	23	UP TO 2	CPOP/XS
SMALL	24.5	UP TO 3.5	CPOP/S
MEDIUM	26	UP TO 5.5	CPOP/M
LARGE	27.5	UP TO 7.5	CPOP/L
X-LARGE	29.5	UP TO 10	CPOP/XL

FOOTSHIELD HEX INSOLE



			SHOE/BOOT	WEDGE
SIZE	SOLE LENGTH CMS	UK SHOE SIZE	PRODUCT CODE	PRODUCT CODE
X-SMALL	24	1.5 - 3	FSHIWS/XSL	FSHIBSB/XSR
SMALL	25.7	3.5 - 5	FSHIWS/SL	FSHIBSB/SR
MEDIUM	27.3	5.5 - 7	FSHIWS/ML	FSHIBSB/MR
LARGE	29	7.5 - 9	FSHIWS/LL	FSHIBSB/LR
X-LARGE	30.7	9.5 - 11	FSHIWSXLL	FSHIBSB/XLR

FOOTSHIELD II SHOE

FEATURES & BENEFITS

- Evidence based rocker sole*
- 10 degree toe out angle
- Apex of rocker at 60% of shoe length
- 20 degree rocker angle
- Removable memory foam insole
- Square-toe design to accommodate the presence of K wires
- Touch-and-close closure, helps to minimise pressure on the forefoot, yet is adjustable to accommodate for bulky dressings
- Can be used on its own or with the Hex plantar foot offloading insole
- Removable toe cap
- Available left and right



Optional Hex insole

			WITH HEX INSOLE			
		SOLE LENGTH	LEFT	RIGHT	LEFT	RIGHT
SIZE	CMS	UK SHOE SIZE	PRODUCT CODE	PRODUCT CODE	PRODUCT CODE	PRODUCT CODE
X-SMALL	24	1.5 - 3	FSS-II/XSL	FSS-II/XSR	FSS-IIH/XSL	FSS-IIH/XSR
SMALL	25.7	3.5 - 5	FSS-II/SL	FSS-II/SR	FSS-IIH/SL	FSS-IIH/SR
MEDIUM	27.3	5.5 - 7	FSS-II/ML	FSS-II/MR	FSS-IIH/ML	FSS-IIH/MR
LARGE	29	7.5 - 9	FSS-II/LL	FSS-II/LR	FSS-IIH/LL	FSS-IIH/LR
X-LARGE	30.7	9.5 - 11	FSS-II/XLL	FSS-II/XLR	FSS-IIH/XLL	FSS-IIH/XLR

*Effect of rocker shoe design features on forefoot plantar pressures in people with and without diabetes. Chapman et Al, Clinical Biomechanics. 28(6):679-85, 2013 July



CONTACT US

M.E.D. SURGICAL CUSTOMER SERVICE

Sisk Healthcare UC, t/a M.E.D. Surgical
Howth Junction Business Centre
Dublin 5, D05 HP92