



Rüsch
TruLite Secure
 Single-Use Laryngoscope Blade and Handle Combination

Ideal for OR Intubations, Crash Carts and Emergency Situations

A Complete Solution to Reduce the Risk of Cross-Contamination

Concerns about proper cleaning of tools used for tracheal intubation have been growing as bacteria are becoming more resistant and no perfect decontamination procedures exist.¹ While most concerns about laryngoscope cleaning focus on the blade because of the patient contact, handles are becoming more of a concern as studies have shown that handles also carry bacteria.²

The Rüsch TruLite Secure Single-Use Laryngoscope Blade and Handle Combination from Teleflex is a single-use solution that consists of a combined laryngoscope blade and handle. This innovative system offers the opportunity to reduce the potential for infection associated with reusable devices.

Dependable Design and Construction

The blade is constructed of stainless steel with rounded edges and a bevelled tip, while the handle combines

high-grade, reinforced plastic with an ergonomic grip. The durable construction of these critical components is designed to withstand the forces commonly associated with intubation.

Confident and Convenient Use

The Rüsch TruLite Secure Single-Use Laryngoscope Blade and Handle Combination is a one-piece blade and handle system that minimises the risk of incompatibility while featuring a wide range of colour-coded blades that help to enable clinicians to quickly identify and select the correct sizes for their patients. The Rüsch TruLite Secure Single-Use Laryngoscope Blade and Handle Combination also improves visibility for the clinician by providing strong, focused and reliable illumination from the integrated LED lighting technology, which is powered by pre-loaded batteries in the handle.

Benefits



The Clinician
 Design helps inspire clinician confidence during critical procedures



Your Institution
 Cost-neutral* way to combat intubation infection risk




The Patient
 Promotes positive patient outcomes by reducing the risk of cross-contamination

*Individual costs and savings may vary based on a number of factors, including specific markets, services, and rates.

Rüsch

TruLite Secure

Single-Use Laryngoscope Blade and Handle Combination



Rüsch TruLite Secure Single-Use Laryngoscope Blade and Handle Combination: Fast Facts

Fully disposable solution to help reduce the risk of cross-contamination and eliminate reprocessing costs.

Compact and lightweight. Suitable for use in the OR, crash carts, airway kits and EMS bags.

Reinforced critical components are designed to perform in the toughest OR and field conditions.

2 AAA batteries are pre-loaded in the handle for immediate use in emergency situations. Batteries are validated for a 3-year shelf life.

Bright blue-white LED light source to enhance tissue differentiation.

Rüsch TruLite Secure Single-Use Laryngoscope Blade and Handle Combination

REF.	DESCRIPTION	OVERALL LENGTH (MM)	DISTAL WIDTH (MM)	WEIGHT (WITH BATTERIES) (G)	COLOUR CODE	CASE QTY
MACINTOSH BLADES						
004671002	Mac 2	100	11.3	108	Turquoise	10
004671003	Mac 3	130	15	125	Yellow	10
004671004	Mac 4	155	15	135	Pink	10
MILLER BLADES						
004670010	Mill 00	67	11.8	97	Beige	10
004670000	Mill 0	78	11.8	98	Purple	10
004670001	Mill 1	103	11.8	104	Orange	10
004670002	Mill 2	155	13.7	125	Grey	10
004670003	Mill 3	195	13.7	137	Green	10
004670004	Mill 4	205	13.7	141	Lavender	10

References:

- Galinski M, Adnet F, Tran D, et al. *Eur J Anaesthesiol.* 2003;20(9):731-735.
- Williams D, Dingley J, Jones C, Berry N. *Journal of Hospital Infection.* 2010;74(2):123-128.

Rx only

Contents of unopened, undamaged package are ready for use and for single-patient use. Refer to package insert for current warnings, indications, contraindications, precautions, and Instructions for Use.

Teleflex, the Teleflex logo, Rüsch, and TruLite are trademarks or registered trademarks of Teleflex Incorporated or its affiliates, in the U.S. and/or other countries. Information in this document is not a substitute for the product Instructions for Use. The products in this document may not be available in all countries. Please contact your local representative. All data current at time of printing (07/2017). Subject to technical changes without further notice. © 2017 Teleflex Incorporated. All rights reserved. 94 08 64 - 00 00 01 - REV A - 07 17 PDF



RuschReliability.com



Distributed by:

Teleflex Headquarters International, Ireland · Teleflex Medical Europe Ltd. · IDA Business & Technology Park
Dublin Road · Athlone · Co Westmeath · Tel. +353 (0)9 06 46 08 00 · Fax +353 (0)14 37 07 73 · orders.intl@teleflex.com
United Kingdom Tel. +44 (0)14 94 53 27 61 · info.uk@teleflex.com
South Africa Tel. +27 (0)11 807 4887 · assist.africa@teleflex.com